


LongerPump®硅胶管



管号	内径 (mm)	外径 (mm)	壁厚 (mm)	包装 (米/包)	硬度 (邵A)	扯断强度 (Mpa)	扯断延长 率 (%)	撕裂强度 (kN/m)	材料
13#	0.8	4.2	1.7	2/5/15	50-55	≥7.0	≥300	≥30	铂金硫化硅胶管
14#	1.6	5	1.7	2/5/15	50-55	≥7.0	≥500	≥30	
19#	2.4	5.8	1.7	2/5/15	50-55	≥7.0	≥500	≥30	
16#	3.2	6.4	1.6	2/5/15	50-55	≥7.0	≥500	≥30	
25#	4.8	8	1.6	2/5/15	50-55	≥7.0	≥500	≥30	
17#	6.4	9.6	1.6	2/5/15	50-55	≥7.0	≥500	≥30	
18#	7.9	11.1	1.6	2/5/15	50-55	≥7.0	≥500	≥30	
15#	4.8	9.6	2.4	2/5/15	50-55	≥7.0	≥500	≥30	
24#	6.4	11.2	2.4	2/5/15	50-55	≥7.0	≥500	≥30	
35#	7.9	12.7	2.4	2/5/15	56-60	≥7.0	≥300	≥30	
36#	9.5	14.3	2.4	2/5/15	56-60	≥7.0	≥300	≥30	
73#	9.5	16.1	3.3	2/5/15	50-55	≥7.0	≥500	≥30	
82#	12.7	19.3	3.3	2/5/15	50-55	≥7.0	≥500	≥30	
88#	12.7	22.3	4.8	2/5/15	56-60	≥7.0	≥300	≥30	
92#	25.4	35	4.8	2/5/15	50-55	≥7.0	≥500	≥30	

* 使用温度范围:-30°C~250°C

* 包装标准:单层包装